

1. Download the latest version of ITRANS: <http://www.goitrans.com/itrans-support-downloads/>
2. Prior to installing you will need to see what the **Connection ID** will be.
3. Open ABELDent/File/Update/CDANet/Network Information/ICA-Itrans/Connection ID = **ICA or ITRANS**

The screenshot shows the 'Network Information Dialog' window. On the left, a list of providers is shown with 'ICA - ITRANS' selected. The right side of the dialog has the following fields:

- Network ID: ICA
- Name: ITRANS
- Contact Name: (empty)
- Phone: (empty) Ext: (empty)
- Connection IDs:
  - Sending Transactions: ITRANS
  - Payment Reconciliation: ITRANS
  - Summary Reconciliation: ITRANS
- Outstanding Transactions Support:
  - Center directed?:
  - Network directed?:
- Reconciliation Support:
  - Center directed?:
  - Network directed?:
- Controller Type: CDA
- Network Transaction Version: 04
- Gateway Time Out: 120 sec

Buttons at the bottom include 'Add', 'Delete', 'Ok', and 'Cancel'.

The screenshot shows the 'Network Information Dialog' window. On the left, a list of providers is shown with 'ICA - ITRANS' selected. The right side of the dialog has the following fields:

- Network ID: ICA
- Name: ITRANS
- Contact Name: (empty)
- Phone: (empty) Ext: (empty)
- Connection IDs:
  - Sending Transactions: ICA
  - Payment Reconciliation: ICA
  - Summary Reconciliation: ICA
- Outstanding Transactions Support:
  - Center directed?:
  - Network directed?:
- Reconciliation Support:
  - Center directed?:
  - Network directed?:
- Controller Type: CDA
- Network Transaction Version: 04
- Gateway Time Out: 120 sec

Buttons at the bottom include 'Add', 'Delete', 'Ok', and 'Cancel'.

4. Browse to **C:\ABELDent\Gateway** Create a folder from Step 3 called **(ICA or ITRANS)**
5. Install ITRANS, in the new folder location **C:\ABELDent\Gateway\ICA or ITRANS**
6. After call ITRANS you will need to download a certificate for each provider. 1-866-788-1212

For more information about ITRANS please visit: <http://www.goitrans.com/>