
Moving ABELDent DCG to NEW SERVER

LOG IN AS "ADMINISTRATOR" FOR TASKS BELOW, NOT AN ACCOUNT WITH ADMIN PRIVILEGES

- 1) From **OLD SERVER** copy the "**C:ABELDent Folder**". Place it on the **NEW SERVER**
- 2) Install the **exact same DCG Version of ABELDent** on top the "**ABELDent Folder**". You can update to a newer version of "**DCG ABELDent**" after.
- 3) Open the **ABELDent** program to validate the Data looks correct.
- 6) Reinstall client machines select client **radio button** and indicating new **SERVER NAME**.
- 4) If office uses **Itrans** reinstall using the guide "**How to Install Itrans**"
- 5) Once installed - Contact Itrans and install a new **Digital Certificates**, then reboot the server.
- 6) Call ABELDent support to **re-license** SERVER permanent license @ **1-800-263-5104** ext. **2**